

Ziglite Smart Health Care Products Co., Ltd. 節 亮 康 護 產 品 有 限 公 司



Air Powered 氣動式康復訓練器

Pneumatic Rehabilitation Trainer



CATALOG

Ziglite Smart Health Care Products Co., Ltd. 節 亮 康 護 產 品 有 限 公 司 Your Smart Health Care Solution Partners



www.ziglite.com



Office +852 3107 7500

Fax +852 3544 0462

Whatsapp +852 6827 2010

Email info@ziglite.com

Web www.ziglite.com

Office Hour

Mon - Fri : 9:30 - 5:30pm

Hong Kong Office Addres

Room B3, 18/F Bonsun Industrial Building, 364-366 Sha Tsui Road, Tsuen Wan, Hong Kong



ABOUT US

關於我們

Ziglite Smart Health Care Products Co., Ltd. Incorporated in 2018, we work with our partners to bring their products to the local market. Mainly engaged in the management of health care products, the main customers are hospitals and clinics, homes for the elderly and home patients.

We are committed to the development of excellent health care products, dedicated to the needs of a wide range of different customers, get the best care in the family treatment process, and greatly improve the effect of rehabilitation. Get the best care during the family treatment process and greatly improve the healing effect. We will provide quality products with stable quality and stable performance for customers.

We specialize in a wide range of products including: anti-virus and air freshening, incontinence care, mobility aids, decompression supplies, etc.

公司於2018年成立,我們與合作夥伴一起把他們的產品推向本地市場。以經營康護用品為主,主要客戶為醫院及診所,安老院及家居病人等。

我們致力於發展卓越的康護醫療產品,竭誠照顧廣大不同客戶的需要,在家庭治療過程中得到最好的護理,大大提高康復治療效果。在家庭治療過程中得到最好的護理,大大提高康復治療效果。務求為客戶提供品質優越及性能穩定的醫療產品。

我們專營的產品種類繁多,包括:抗病毒及空氣清新,失禁護理,助行用品,減壓用品等。





目录 Content

氣動訓練功能&原理 Function & Principle Of Pneumatic Therapy

氣動式康復訓練 Pneumatic rehabilitation ther	ару	P2
氣動類產品 Pneumatic products	•••••	P4
空氣波壓力按摩 Air wave pressure therapy	•••••	P6
空氣波類產品 Air wave products		P8



氣動康複訓練系列

Pneumatic Rehabilitation Training Series



Host + Air Valve Training Gloves (Lite Version)

主機+ 氣閥訓練手套 精簡款



Host + Air Valve Training Gloves (Upgraded)

主機+氣閥訓練手套+鏡像手套 升級款



Host + Air Valve Training Gloves + Mirror Gloves (Hot Compress & Pulse)

主機+氣閥訓練手套+鏡像手套 熱敷&脈衝款



Host + Training Gloves + Mirror Gloves + External Control Handle

(Professional Edition)

主機+訓練手套+鏡像手套+外控手柄 專業版



Host + Wrist Training Gloves (Wrist Training Gloves)

主機+分指腕關節訓練器 腕關節訓練



Host + Training Gloves + Mirror Gloves + Massage set

(comprehensive Training Package) 主機+訓練手套+鏡像手套 +按摩5件套 綜合訓練組合套



Host + Wrist Training Gloves (Wrist Training Gloves)

主機+分指腕關節訓練器 腕關節訓練



空氣波按摩系列 **Air Wave Massage Product Series**



Host + Six Chamber Massage Coat For Double Upper Limb

主機+雙隻 六腔手臂按摩套



Host + Six Chamber Massage Coat For Double Lower Limb

主機+雙隻 六腔腿部按摩套



Host + Six Chamber Massage Coat For Massage

主機+六腔腰部按摩套



Host + Four Chamber Massage Coat For Single Upper Limb

主機+單隻四腔手臂按摩套



Host + Four Chamber Massage Coat For Single Lower Limb

主機+單隻四腔腿部按摩套



Host + Four Chamber Waist Massage Coat

主機+四腔腰部按摩套



Controller + Thigh Massage Cover + Calf Foot Cover 控制器+大腿按摩套+小腿連腳套



主機+腿部按摩套



Principle Of Pneumatic Rehabilitation Tra

理

02

PART

氣動康復功能實現&原理

Function & Principle Of Pneumatic Rehabilitation Training

Pneumatic 氣動式康復訓練

Pneumatic Rehabilitation Training

氣動康復訓練原理 Principles of Pneumatic Rehabilitation Training

氣動式關節智能康復系統,是以空氣為動力,自動驅使 手指手腕等關節被動訓練,可有效降低肌張力、改善血液迴 圈、緩解關節水腫和增加關節活動度。

又能通過訓練和運動想像等反射腦部神經,促進腦部神經康復。是目前針對手指痙攣、麻痹、癱瘓等症狀有效的康復訓練設備。

Pneumatic joint intelligent rehabilitation system, which is driven by air and automatically drives the passive training of joints such as fingers and wrists, can effectively reduce muscle tension, improve blood circulation, alleviate joint edema and increase joint mobility.

It can also reflect brain nerves through training and motor imagination to promote brain nerve rehabilitation. It is currently an effective rehabilitation training device for symptoms such as finger spasms, paralysis, and paralysis.



康復功能實現 Realization Of Rehabilitation Function

康復是一項漫長又複雜的過程,我們推出的氣動式康 復訓練儀目的是作為補充手指恢復過程中的一個重要階 段的康復訓練器具。

利用短促等長練習TENS法則的訓練模式讓患者體會手指肌肉被動訓練的過程,並進行初級主動訓練,更早的進入下一個階段訓練。

根據患者疾患情況,通過設定調節康復儀指套活動壓力,被動地驅使手指/手腕/腳踝活動,不但可使手指/手腕/腳踝痙攣、麻痹、癱瘓等症狀得到極大改善,而且還可以間接作用於腦部神經及血管,來促進腦部損傷或神經損傷術後的康復。

Rehabilitation is a long and complex process. Our pneumatic rehabilitation training instrument aims to supplement the rehabilitation training instrument at an important stage in the process of finger recovery.

Using the training mode of short equal length exercise tens rule, let the patients experience the process of passive training of finger muscles, carry out primary active training, and enter the next stage of training earlier.

According to the patient's disease situation, by setting and adjusting the finger cuff activity pressure of the rehabilitation instrument and passively driving the finger / wrist / ankle activity, it can not only greatly improve the symptoms of finger / wrist / ankle spasm, paralysis and paralysis, but also indirectly act on the brain nerve and blood pipe to promote the rehabilitation after brain injury or nerve injury.


<u>Q</u> **Rehabilitation Training**

產品特色

Product Features

- 以空氣壓力作為動力,自動驅使手指/手腕被動訓 練、操作簡單、安全。
- 手指、手腕、腳腕等可作用於多個部位一機多用高性
- 被動訓練、主從式鏡像訓練結合,以達到手指功能康 復和腦功能重建雙重目的。
- 伸展位持續時間可調,壓力可調,以適用不同肌張力 患者。
- 可以促進上肢和手的功能,強化患者的自信心,輔助心 裏治療。
- 設備設計精巧,可在床邊及任意空間使用。
- 設備操作簡便,家庭人員從旁協助即可,節省醫護成 本。

- With air pressure as power, it automatically drives fingers / wrists for passive training, which is simple and safe.
- Fingers, wrists, ankles, etc. can act on multiple parts, one machine with multi-purpose and high
- Passive training and master-slave mirror training are combined to achieve the dual purpose of finger function rehabilitation and brain function reconstruction
- The duration of the stretch position is adjustable, and the pressure is adjustable to suit patients with different muscle tensions.
- It can promote the function of upper limbs and hands, strengthen patients' self-confidence and assist psychological treatment.
- The equipment is exquisitely designed and can be used beside the bed and in any space.
- The equipment is easy to operate, and family personnel can assist from the side to save medical





pneumatic rehabilitation training

適用範圍 Scope Of Application

- 腦卒中、腦癱、腦血栓、脊髓損傷等中樞神經系統損 傷引起的肢體功能障礙。
- 骨科手術、外傷等外周神經損傷引起的肢體功能障
- 手指功能障礙從軟癱期至恢復期的康復需求的患 者。
- Limb dysfunction caused by central nervous system injury such as stroke, cerebral palsy, cerebral thrombosis and spinal cord injury.
- Limb dysfunction caused by peripheral nerve injury such as orthopedic surgery and trauma.
- Patients with finger dysfunction who need rehabilitation from the stage of soft paralysis to the recovery stage.

V SKD ZIGLITE ▲ SKB

通過氣泵控制訓練終端進行康復訓練運動。

其原理為:主機氣泵將氣壓連接到訓練手套的氣囊 內,並控制氣泵打氣和抽氣工作來控制訓練終端的氣囊的 擴張和收縮從而實現手指、手腕、腳腕的來回運動。

電源為寬電壓設計,相容AC100-240V/50-60Hz,全 球通用。顯示屏幕使用LED屏,顯示內容更清晰自然。

操作介面簡單明白,更易學易操作,帶中/英文語言版 本選擇,適用性更強。

Training Gloves



手指訓練 / Finger training

以空氣壓力作為動力,通過主機控制自動驅使手 指抓、握、伸展等被動訓練。

有效的抑制手指的痙攣,增加手指關節活動度, 改善血液迴圈,達到放鬆手指關節的效果。

With the air pressure as the power, the passive training such as finger grasping, grasping and stretching is automatically driven by the host control.

It can effectively inhibit finger spasm, increase finger joint activity, improve blood circulation and relax finger joints.

Products

氣動康復訓練產品

Pneumatic Rehabilitation Training

康復訓練器主機

Rehabilitation Device Host

實

現

原

理

04

Principle


Qf

T

neum

atic

Rehabilitation



The rehabilitation training is carried out through the air pump control training terminal.

The principle is: the host air pump connects the air pressure to the air bag of the training glove, and controls the pumping and pumping of the air pump to control the expansion and contraction of the air bag of the training terminal, so as to realize the back and forth movement of fingers, wrists and ankles.

The power supply is designed for wide voltage, compatible with AC100-240V/50-60Hz, universally applicable. The display screen uses LED screen, the display content is more clear and natural.

The operation interface is simple and clear, easier to learn and operate, with Chinese / English language version selection and stronger applicability.





氣閥訓練手套:氣閥單獨控制每根手指 Training gloves: each finger is individually controlled by an air valve switch

實

現

理

06



▲ SA-05 訓練手套:多根氣管控制, 搭配主機可進行分指控制

Training gloves: multiple trachea control, It can be controlled separately with the host



鏡像手套、鏡像手柄,控制鏡像訓練 Mirror gloves, mirror handles, control mirror training



手腕訓練 / Wrist training

以空氣壓力作為動力,通過主機控制自動驅使手腕活動被動訓練。幫助手腕往上抬起,可幫助腕關節受阻的 患者進行拉伸,緩解腕部筋膜的僵硬

With the air pressure as the power, the wrist is automatically driven by the host control for passive training.

Helping the wrist lift up can help patients with wrist joint obstruction stretch and alleviate the stiffness of wrist fascia



▲ SA-06

分指手腕關節訓練器 Split finger wrist joint trainer

手腕訓練手套(基礎款) Wrist Training Gloves (Basic)

手&腳部氣壓按摩

Hand & foot air wave massage

以空氣壓力作為動力,控制手/腳套的氣囊充氣和抽氣狀態來控制手/腳套內的氣囊的擴張和收縮,從而降低肌張力、改善血液循環、緩解水腫。

Take the air pressure as the power to control the inflation and pumping state of the air bag in the hand / foot cover to control the expansion and contraction of the air bag in the hand / foot cover, so as to effectively reduce muscle tension, improve blood circulation and alleviate edema.



▼ SA-01 空氣波按摩手套 Air wave massage gloves



⋖ SA-02

空氣波足底按摩腳套套 Air wave massage foot cover

Air Wave 空氣波壓力按摩

Air Wave Pressure Massager

空氣波壓力按摩原理 Principle Of Air Wave Pressure Massage

空氣波壓力按摩又稱循環壓力按摩、四肢循環儀或壓力抗 栓泵,物理療法。

空氣波壓力按摩儀主要通過對多腔氣囊有順序的反復充放氣,形成了對肢體和組織的循環壓力,對肢體的遠端到肢體的近端進行均勻有序的擠壓,促進血液和淋巴的流動,改善微循環的。加速肢體和組織血液回流,有助於預防血栓的形成、預防肢體水腫,能夠直接或間接治療與血液淋巴循環相關的諸多疾病。

Air wave pressure massager, also known as circulatory pressure massage, limb circulatory device or pressure anticoagulant pump, physical therapy.

The air wave pressure massager instrument mainly forms the circulating pressure on the limbs and tissues through the sequential repeated inflation and deflation of the multi cavity air bag, and evenly and orderly compresses the distal end of the limbs to the proximal end of the limbs, so as to promote the flow of blood and lymph and improve the microcirculation. Accelerating the blood reflux of limbs and tissues is helpful to prevent thrombosis and limb edema. It can directly or indirectly treat many diseases related to blood lymphatic circulation.



康復功能實現 Realization Of Rehabilitation Function

通過被動均勻的按摩作用,隨著血液循環的加速。可以加速血液中代謝廢棄物,炎症因數、致痛因數的吸收。可以防止肌肉萎縮,防止肌肉纖維化,加强肢體的含氧量,有助於輔助解决因血液循環障礙引起的疾病(如骨股頭壞死等)。

Through passive and uniform massage, with the acceleration of blood circulation. It can accelerate the absorption of metabolic wastes, inflammatory factors and pain causing factors in the blood. It can prevent muscle atrophy, prevent muscle fibrosis, strengthen the oxygen content of limbs, and help to solve the diseases caused by blood circulation disorders (such as ring death of femoral head, etc.).

Principle

Qf

Pneumatic

Rehabilitation

Training

現

理

08


Training


產品特色

Product Features

- 以空氣壓力作為動力,設備操作越來越簡單,醫用家用 皆可以。
- 安全、綠色、無創傷,符合現代醫學發展方向。
- 氣壓按摩,放鬆舒適。
- 對某些疾病具有多重功效。
- 簡便,無需護理人員從旁協助節省醫護成本。

- With air pressure as power, the operation of treatment equipment is becoming simpler and simpler, which can be used for medical and household purposes.
- Safe, green and non-invasive, in line with the development direction of modern medicine.
- Pneumatic massage, relaxation and comfort.
- It has multiple effects on some diseases.
- Simple and convenient, without the assistance of nursing staff, saving medical costs.





適用範圍 Scope Of Application

● ト、下肢體水腫:

上、下肢體的原發性和繼發性淋巴水腫,慢性靜脈源性 水腫,脂性水腫,混合性水腫等。

治療原理是促進血液循環,淋巴循環,將一些引起疼痛, 不舒服的代謝產物和炎性致痛物質都擠到主循環裏清 除掉,從而消除水腫。

偏癱,截癱及癱瘓的病人:

偏癱、截癱、癱瘓、長期臥床的患者,因血流緩慢,肌肉無 收縮力。應用空氣波壓力治療儀,反復對肢體進行加壓 後再卸壓,從而產生如同肌肉的收縮和舒張作用,促進 靜脈血液和淋巴液循環,並取得充分按摩效果。預防深 靜脈血栓和防治下肢肌肉萎縮。

糖尿病足,糖尿病末梢神經炎:

增加了神經血流灌注和氧合作用,並且提高神經的氧耗 量從而達到改善功能的目的,减轻肢体麻木症狀。

肢體血液循環不好,靜脈功能不全,中老年人:

促進血液循環,解除肌肉疲勞,緩解疼痛,麻痹不適肢體 的康復。

空氣波壓力按摩儀運用間歇壓力,通過空氣波的反復膨 脹和收縮作用,可以達到改善血液循環,擴展活化血管 的效果。有助於抗血栓形成和改善循環,清除血液中代 謝廢棄物,加強肢體氧合度。

Edema of upper and lower limbs:

Primary and secondary lymphedema, chronic venous edema, fatty edema, mixed edema, etc. of upper and lower limbs.

Promote the metabolism of inflammatory pain causing substances in the main circulation, so as to eliminate edema.

Patients with hemiplegia, paraplegia and

The air wave pressure therapeutic instrument is used to repeatedly pressurize the limbs and then relieve the pressure, so as to produce the contraction and relaxation of muscles, prevent deep vein thrombosis and lower limb muscle atrophy.

Diabetic foot, diabetic peripheral neuritis:

It increases the nerve blood perfusion and oxygen cooperation, and increases the oxygen consumption of the nerve, so as to improve the function and reduce the symptoms of limb numbness.

Poor limb blood circulation, venous insufficiency, middle-aged and elderly people:

Promote blood circulation, relieve muscle fatigue, relieve pain, and recover paralyzed and uncomfortable limbs


It helps to prevent thrombosis and improve circulation, remove metabolic waste from blood and strengthen limb oxygenation.

Products 空氣波壓力按摩產品

Air Wave Pressure Massage

康復治療器主機

Therapeutic Device Host



主機氣泵將氣壓連接到肢體氣囊內,控制氣泵的加 壓和泄壓工作狀態來控制多腔氣囊有順序的反復充放氣, 形成了對肢體和組織的循環壓力,促進血液和淋巴的流 動,改善微循環。

The air pump of the host connects the air pressure to the limb airbag, controls the pressurization and pressure relief working state of the air pump to expand and contract the airbag of the upper limb sleeve, forms the circulating pressure on the limbs and tissues, promotes the flow of blood and lymph and improves microcirculation.

▲ M-01

Airbag Massage Terminal



上肢按摩 / Upper limb massage

應用空氣波壓力按摩儀,反復對肢體進行加壓後再卸壓,促進血液 循環,淋巴循環,緩解手臂肌肉疲勞,减少手臂僵硬等問題。

Apply air wave pressure massager instrument to repeatedly pressurize the limbs and then relieve the pressure, promote blood circulation and lymphatic circulation, alleviate arm muscle fatigue and reduce arm stiffness.



六腔氣動手臂按摩套 Pneumatic arm massage coat - six chamber





四腔氣動手臂按摩套 Pneumatic arm massage coat - four chambers



腰部按摩 / Wrist massage

應用空氣波壓力按摩儀,反復對腰部進行加壓後再卸壓,促進血液循環,淋巴循環,緩解腰部疲勞酸楚疼痛等問題。

he air wave pressure massager is used to repeatedly pressurize the waist and then relieve the pressure to promote blood circulation, lymphatic circulation, and relieve waist fatigue and pain.



▲ M-65


六腔氣動腰部按摩套 Pneumatic waist massage cover - six chambers





四腔氣動腰部按摩套 Pneumatic waist massage cover

⋖ M-45





下肢按摩 / Lower limb massage

應用空氣波壓力按摩儀,反復對下肢進行加壓後再卸壓,促進血液循環,淋巴循環,改善腿麻、脹、腫、疼等問題,輔助術後恢復,預防靜脈曲張、血栓等問題。

Apply air wave pressure massager instrument to repeatedly pressurize the limbs and then relieve the pressure, promote blood circulation and lymphatic circulation, improve leg numbness, swelling, swelling and pain, assist postoperative recovery, and prevent varicose veins, thrombosis and other problems.





空氣波壓力腿部按摩儀 Air wave pressure leg massage instrument



空氣波壓力腿部按摩儀 Air wave pressure leg massage instrument



複

實

現

原

理

10

Principle Of Pneumatic

Rehabilitation Training

Pneumatic Training Product Series

氣動康複訓練系列

Pneumatic Training Product Series

手指伸展訓練器(精簡版)

Finger Stretching Trainer (Lite Version)

Model: SKA

主機+氣閥訓練手套

Host + Air Valve Training Gloves

▶ 產品特點/Features

- 全自動化被動訓練模式,促進上肢和手的功能,並能可通過單個 手指氣閥對單個或指定的手指進行定制康復訓練功能。
- 使用軟矽膠+彈力面料,訓練過程中不會造成患肢的二次傷害。
- 訓練手套的每根手套氣管上都帶單個相互獨立的氣閥開關,可 通過該開關的打開或關閉來進行手指針對性選擇訓練。
- 訓練力量可設置,滿足不同患者的需求。
- 可居家訓練,保持醫療機構手康復訓練成果,並能促進手指康復
- The fully automatic passive training mode promotes the functions of upper limbs and hands, and can customize the rehabilitation training function of single or designated fingers through a single finger air valve.
- The use of soft silicone rubber + elastic fabric will not cause secondary injury to the affected limb during training.
- Each glove can be trained by opening and closing the air valve on each finger independently.
- The training strength can be set to meet the needs of different
- It can train at home, maintain the results of hand rehabilitation training in medical institutions, and promote the finger rehabilitation journey.



▶ 產品規格/Specifications

主機 / Host

材質 / Material 鋁合金外殼 aluminum alloy shell

電源 / Power Supply input AC100-240V, output DC12V, 50/60Hz

功率 / Power max 12W, min 6W, standby 0.6W

訓練模式 / Training mode 被動訓練,鏡像訓練 passive training, mirror training

設置 / Setting 力度:1~9檔

工作時間 / Time Strength setting: 1~9 15min

運行噪音 / Operating noise ≤60 dB

尺寸 / Size (L*W*H) 重量 / Weight 103*90*60mm 0.47kg

PVC+soft silicone+elastic cloth

左手,右手:S,M,L,XL,XXL號 left and right:S,M,L,XL,XXL

重量 / Weight 0.22kg

氣閥訓練手套 / Air Valve Training Gloves

Training Product Series

列

12

Pneumatic

Serie

手指伸展訓練器(升級款)

Finger Stretching Trainer (Upgraded)

Model: SKB

主機+氣閥訓練手套+鏡像手套

Host + Air Valve Training Gloves + Mirror Gloves

▶ 產品特點/Features

- 全自動化被動訓練模式+鏡像訓練,有助於促進患手肌力恢復, 促使主動屈伸運動出現。
- 高清LED數碼管,數據顯示清晰。
- 所有功能鍵均有圖示顯示,觸摸按鍵操作更輕鬆便捷。
- 操作界面簡單明白,更易學易操作,可居家訓練。
- 帶運行蜂鳴器工作音提示,隨時掌握訓練情況。
- "訓練力量"及"訓練週期"可設定,滿足不同患者的需求。
- The fully automatic passive training mode+mirror training is helpful to promote the recovery of muscle strength of the affected hand and promote the emergence of active flexion and
- High-definition LED digital tube, clear data display.
- All function keys are shown in the figure, which makes the operation of touch keys easier and more convenient.
- The thumb and the other four fingers are controlled separately
- The operation interface is simple and clear, easy to learn and operate, and can be trained at home.
- With the active buzzer working sound prompt, you can master the training situation anytime.
- "Training strength" and "training cycle" can be set to meet the needs of different patients,



▶ 產品規格/Specifications

主機 / Host

材質 / Material 塑膠 plastic

電源 / Power Supply input AC100-240V, output DC12V, 50/60Hz

功率 / Power

max 10W, min 8W, standby 2W

訓練模式 / Training mode

被動訓練,鏡像訓練 passive training, mirror training

設置 / Setting

力度:1~9檔,週期:1~9檔

Strength setting: 1~9 Training cycle: 1~9

運行噪音 / Operating noise ≤60 dB

尺寸 / Size (L*W*H) 251*204*111mm

氣閥訓練手套 / Air Valve Training Gloves

PVC+軟硅膠+透氣彈力面料 PVC+soft silicone+breathable stretch fabric

左手,右手:S,M,L,XL,XXL號 left and right:S,M,L,XL,XXL

鏡像手套 / Mirror Gloves

材質 / Material 棉線 cotton thread

尺碼 / Size

左手,右手:S,M,L 號 left and right:S, M, L



列

14

Pneumatic

ing

Product

S

熱敷低頻脈衝手功能訓練器

Hot Compress & Low-frequency Pulse **Hand Trainer**

Model: SKC

主機+氣閥訓練手套+鏡像手套

Host + Air Valve Training Gloves + Mirror Gloves

▶ 產品特點/Features

- 全自動化被動訓練模式+鏡像訓練,有助於促進患手肌力恢復, 促使主動屈伸運動出現。
- 高清液晶顯示幕,功能清晰標明,老人、病患可獨立操作。
- 新增了低頻脈衝、熱敷加熱兩個實用功能,讓手部訓練更加科學
- 熱敷3檔溫度調節:45°/52°/62°三檔溫度。
- 低頻脈衝錘擊、按摩、混合三種模式可選,通過低頻脈衝直達神 經促進運動,有效緩解肌張力,更好的達到訓練效果。
- "訓練力量"及"訓練週期"可設定,滿足不同患者的需求。
- The fully automatic passive training mode+mirror training is helpful to promote the recovery of muscle strength of the affected hand and promote the emergence of active flexion and extension movement.
- High definition LCD display screen with clear function indication, capable of independent operation for the elderly and patients.
- Added two practical functions: low-frequency pulse and hot compress heating, making hand training more scientific and efficient.
- Temperature adjustment for hot compress in three levels: 45 ° /52 °
- There are three modes to choose from: low-frequency pulse hammering, massage, and mixing. The low-frequency pulse directly reaches the nerves to promote movement, effectively relieving muscle tension and achieving better training results.
- "Training strength" and "training cycle" can be set to meet the needs of different patients.

▶ 產品規格/Specifications

主機 / Host

材質 / Material 電源 / Power Supply 塑膠 plastic input AC100-240V, output DC7.4V, 50/60Hz

功率 / Power max 21W, min 3.5W, standby 1W

訓練模式 / Training mode 被動訓練,鏡像訓練 passive training, mirror training

設置 / Setting 力度:1~9檔,週期:1~9檔

Strength setting: 1~9 Training cycle: 1~9

訓練時間設定 / Training time

20min/30min/40min The default is: 20 min 運行噪音 / Operating noise 尺寸 / Size (L*W*H)

≤60 dB 135*120*60mm

氣閥訓練手套 / Air Valve Training Gloves

PVC+軟硅膠+透氣彈力面料 PVC+soft silicone+breathable stretch fabric

左手,右手:S,M,L,XL,XXL號 left and right:S,M,L,XL,XXL

鏡像手套 / Mirror Gloves

材質 / Material 棉線 cotton thread

尺碼 / Size

左手,右手:S,M,L 號 left and right:S, <u>M,</u> L

品 系

列

16

手指伸展訓練器(專業版)

Individual Finger Stretching Trainer (Professional Edition)

Model: SKD

主機+訓練手套+鏡像手套+外控手柄

Host + Training Gloves + Mirror Gloves + External Control Handle

▶ 產品特點/Features

- 全自動化被動訓練模式+鏡像訓練,有助於促進患手肌力恢復, 促使主動屈伸運動出現。
- 高清LED數碼管,數據顯示清晰。
- 所有功能鍵均有圖示顯示,觸摸按鍵操作更輕鬆便捷。
- 增加了壓力顯示外控手柄,可以清晰顯示訓練中的張開及收縮 的壓力狀況。
- 帶運行蜂鳴器工作音提示,隨時掌握訓練情況。
- "訓練力量"及"訓練週期"可設定,滿足不同患者的需求。
- The fully automatic passive training mode+mirror training is helpful to promote the recovery of muscle strength of the affected hand and promote the emergence of active flexion and extension movement.
- High-definition LED digital tube, clear data display.
- All function keys are shown in the figure, which makes the operation of touch keys easier and more convenient.
- The external control handle for pressure display is added, which can clearly display the pressure of opening and contracting during training.
- With the active buzzer working sound prompt, you can master the training situation anytime.
- "Training strength" and "training cycle" can be set to meet the needs of different patients.



主機 / Host

塑膠 plastic

電源 / Power Supply

功率 / Power

訓練模式 / Training mode

被動訓練,鏡像訓練 passive training, mirror training

設置 / Setting

力度:1~9檔,週期:1~9檔

Strength setting: 1~9 Training cycle: 1~9

≤60 dB

尺寸 / Size (L*W*H)

▶ 產品規格/Specifications

材質 / Material

input AC100-240V, output DC12V, 50/60Hz

max 15W, min 7W, standby 2W

運行噪音 / Operating noise

251*204*111mm

訓練手套 / Training Gloves

材質 / Material PVC+軟硅膠+透氣彈力面料 PVC+soft silicone+breathable stretch fabric

左手,右手:S,M,L,XL,XXL號 left and right: S, M, L,XL,XXL

鏡像手套 / Mirror Gloves

材質 / Material 棉線 cotton thread

尺碼 / Size

左手,右手均碼 left and right,free size

壓力外控手柄

Pressure External Control Handle

材質 / Mate<u>rial</u> 塑膠 Plastic

尺碼 / Size(L*W*H) 140*80*38mm

壓力設置範圍 / Weight

Systolic pressure: 0-100KPa (± 10%) Opening pressure: 0-60KPa (± 10%)

分指手腕關節訓練儀

Split Finger Wrist Joint Training Instrument

Model: SKE

主機+分指手腕關節訓練器

Host + Wrist Joint Trainer

▶ 產品特點/Features

- 採用高分子材料和人體工學設計,舒適耐用
- 拉伸肌肉筋腱,恢復神經運動功能
- 促進局部血液循環,防止肌肉萎縮
- 促使肌肉功能恢復,緩解腕部肌肉僵硬
- 背曲訓練:適用於腕關節術後康復期間關節部位肌肉萎縮,韌帶
- 掌曲訓練:適用於腕關節僵硬彎曲困難患者
- "訓練力量"及"訓練週期"可設定,滿足不同患者的需求
- High polymer material and ergonomic design, comfortable and
- Stretch muscles and tendons to restore neuromotor function.
- Promote local blood circulation and prevent muscle atrophy.
- Promote muscle function recovery and relieve wrist muscle
- Back bending training: applies to the joint's muscular atrophy and ligament laxity during the rehabilitation after wrist surgery.
- Palm flexion training: suitable for patients with wrist stiffness and bending difficulties.
- "Training strength" and "training cycle" can be set to meet the needs of different patients.



▶ 產品規格/Specifications

主機 / Host

材質 / Material 塑膠 plastic

電源 / Power Supply input AC100-240V, output DC12V, 50/60Hz

功率 / Power

max 13W, min 9W, standby 4.5W

訓練模式 / Training mode 被動訓練 passive training

設置 / Setting

力度:1~9檔,週期:1~9檔

Strength setting: 1~9 Training cycle: 1~9

運行噪音 / Operating noise ≤60 dB

尺寸 / Size (L*W*H) 251*204*111mm

重量 / Weight 1.2kg

分指手腕關節訓練器 / Split Finger Wrist Joint Trainer

材質 / Material 高分子材料 polymer materials

尺碼 / Size

左手,右手均碼可調節 left hand, right hand, one size adjustable

重量 / Weight 1.3kg

訓練手套 Air Valve Training Gloves

PVC+軟硅膠+牛奶絲彈力面料 PVC+soft silicone+milk silk stretch fabric

左手,右手:S,M,L號 left and right:S,M,L

重量 / Weight 0.32kg

多功能手部訓練器

Multifunctional Hand Trainer

Model: SKF

主機+訓練手套+鏡像手套+刺猬按摩球

- +空氣波按摩手套+低頻脈衝貼片
- +低頻脈衝手握棒+低頻脈衝穴位鞋



Host + Training Gloves + Mirror Gloves + Hedgehog Massage Ball

- + Air Wave Massage Gloves + Low Frequency Pulse Patch
- + Low Frequency Pulse Hand Grip + Low Frequency Pulse Acupoint Shoes

▶ 產品特點/Features

- 全自動化被動訓練模式+鏡像訓練,通過鏡像運動刺激鏡像神經 元,以促進腦功能及手功能的重塑及恢復。
- 分指訓練模式被動模式和鏡像模式均可使用,可選擇性單指訓 練亦可被動機器帶動手部訓練,單指、循環手指訓練。
- 全觸屏操作,手掌螢幕顯示,單獨選擇需要訓練的手指。
- 訓練手套柔性材質,訓練過程中不會造成患肢的二次傷害。
- The fully automated passive training mode + mirror training stimulates mirror neurons through mirror motion to promote the remodeling and recovery of brain function and hand function.
- The split finger training mode, passive mode and mirror mode can be used. It can be used for selective single finger training, passive machine with hands training, single finger and circular finger training.
- Full touch screen operation, palm screen display, separate selection of fingers to be trained.
- The training gloves are made of flexible material, which will not cause secondary injury to the affected limb during training.

▶ 產品規格/Specifications

主機 / Host

材質 / Material

鋁合金外殼+PVC aluminum alloy shell + PVC

電源 / Power Supply

input AC100-240V, output DC12V, 50/60Hz

功率 / Power

35W

訓練模式 / Training mode

被動單指、對指、多指訓練,鏡像訓練,空氣波按摩,

passive single finger, opposite finger and multi finger training, mirror training, air wave massage, low frequency pulse massage

屈縮和伸張速率 / Flexion and extension rate 1st ~ 5th gear

尺寸 / Size (L*W*H) 230*150*100mm

重量 / Weight 3.2ka

空氣波按摩手套 鏡像手套 / Mirror Gloves

材質 / Material 潛水布料 diving cloth

左手,右手:S,M,L號 left and right:S,M,L

重量 / Weight

Air Wave Massage Gloves

海綿雙面復合邊倫布,氣囊為TPU sponge double-sided composite edge Lunbu, the airbag is TPU

尺碼 / Size(L*W*H) 295*260*105mm

壓力範圍 / Pressure Range 26kap

重量 / Weight



訓

品

系

列

18

Pneumatic

Training

Prod

uct

S

erie



手腕訓練手套 / Wrist Training Gloves

通過主機帶動訓練器的起伏,幫助手腕往上 抬起可幫助腕關節受阻的患者,進行拉伸,緩 解腕部筋膜的僵硬。

The host drives the ups and downs of the trainer to help lift the wrist up, which can help patients with blocked wrist joint stretch and alleviate the stiffness of wrist fascia.

空氣波按摩手套 Air Wave Massage Gloves

應用手部空氣按壓方式,對手部反復按摩治 療,使僵硬手部組織軟化和加快血液流通,再 配合訓練手套治療,療效更快更好。

The application of hand air compression and repeated hand massage can soften the stiff hand tissue and speed up the blood circulation. Combined with training gloves. the curative effect is faster and better.

低頻脈衝穴位鞋 Low Frequency Pulse Acupoint Shoes

通過低頻脈衝對足底穴位進行刺激,類比足 療按摩的管道刺激足底神經,舒緩足底的緊

Stimulate plantar acupoints through lowfrequency pulses, simulate foot massage, stimulate plantar nerves and relieve plantar

手腕訓練手套 / Wrist Training Gloves

材質 / Material ABS塑料,尼龍軟粘帶 ABS plastic, nylon soft adhesive tape

尺碼 / Size(L*W)

左手,右手 left and right:210*290mm

壓力範圍 / Pressure Range 正壓180kap,負壓80kap P+180kap,P-80kap 0.04kg

重量 / Weight 0.38kg

刺猬按摩球 / Hedgehog Massage Ball

材質 / Material

ABS plastic, nylon soft adhesive tape

尺碼 / Size Ø70mm

重量 / Weight

低頻脈衝手握棒 Low Frequency Pulse Hand Grip 通過凸點低頻脈衝刺激對僵硬的手部

神經元進行舒緩按摩,緩解手部筋膜的 緊張感,從而緩解僵硬。

The stiff hand neurons are soothed and massaged by bump low-frequency pulse stimulation to relieve the tension of the hand fascia, so as to alleviate the stiffness.



低頻脈衝貼片 Low Frequency Pulse Patch

將貼片貼在相應的穴位或部位,主機產 生的低頻脈衝會通過貼片進入人體,針 對肌肉組織進行放鬆,緩解筋膜的僵 硬,對全身關節痛點起到按摩舒緩的作

Stick the patch on the corresponding acupoints or parts, and the lowfrequency pulse generated by the host will enter the human body through the patch to relax the muscle tissue, relieve the stiffness of fascia, and massage and relieve the pain points of joints of the whole body.

低頻脈衝貼片 Low Frequency Pulse Patch

材質 / Material 導電凝膠 conductive gel

尺碼 / Size(L*W) 52*52mm*2pcs

重量 / Weight 0.005kg*2pcs

低頻脈衝手握棒

Low Frequency Pulse Hand Grip

導電矽膠 conductive silicone adhesive

尺碼 / Size(W*H*L) Ø35*85mm*2pcs

重量 / Weight 0.06kg*2pcs

低頻脈衝穴位鞋

Low Frequency Pulse Acupoint Shoes

材質 / Material plastic + magnet

尺碼 / Size(L*W*H) M: 255*115*10mm x 2pcs

重量 / Weight

0.15kg*2pcs

分指手腕關節訓練儀

Split Finger Wrist Joint Training Instrument

Model: SKG

主機+分指腕關節訓練器

Host + Wrist Joint Trainer

▶ 產品特點/Features

- 底部彈簧可移動臂托設計,可以隨著鍛煉的過程而移動,方便靈
- 中空設計,符合人體工程學設計,在掌屈/背屈過程中防止擠壓
- 腕管,保護腕管。
- 手指手掌手腕部位均採用開放式的魔術貼設計, 鬆緊可調, 貼合 手腕,分指固定,不分左右,適合更多人手型。
- 四角吸盤,黏貼牢固。
- The bottom spring movable arm support design can move with the exercise process, making it convenient and flexible without
- Hollow design, ergonomically designed to prevent compression of the wrist tube during palm/back flexion, protecting the wrist
- The fingers, palms, and wrists are all designed with open Velcro straps, with adjustable elasticity and fit to the wrist. They are fixed by dividing the fingers, regardless of left or right, suitable for more hand shapes.
- Four corner suction cups, firmly attached.



▶ 產品規格/Specifications

主機 / Host

材質 / Material 塑膠 plastic

電源 / Power Supply input AC100-240V, output DC12V, 50/60Hz

功率 / Power

max 10W, min 8W, standby 2W

訓練模式 / Training mode

被動訓練,鏡像訓練 passive training, mirror training

設置 / Setting

力度:1~9檔,週期:1~9檔

Strength setting: 1~9 Training cycle: 1~9

運行噪音 / Operating noise

≤60 dB

尺寸 / Size (L*W*H) 251*204*111mm

分指手腕關節訓練器 / Split Finger Wrist Joint Trainer

材質 / Material polymer materials 運行角度 / Operating angle

尺碼 / Size (L*W*H) 340*210*70mm

重量 / Weight 1.315kg



SA-03

SA-05

訓練手套

單管氣閥訓練手套

左手,右手:S,M,L,XL,XXL號

left and right: S, M, L,XL,XXL Single pipe air valve training gloves

配件選配 Optional accessories

材質 / Material

the airbag is TPU

材質 / Material

材質 / Material

尺碼 / Size

PVC+軟硅膠+诱氣彈力面料

左手,右手:S,M,L,XL,XXL號

left and right : S, M, L,XL,XXL

PVC+soft silicone+breathable stretch fabric

26kap

海綿雙面復合邊倫布,氣囊為TPU

尺碼 / Size 295*260*105mm

壓力範圍 / Pressure Range

PVC+軟硅膠+诱氣彈力面料

sponge double-sided composite edge Lunbu




SA-02

空氣波按摩腳套

Foot cover

Air wave massage

Air valve training gloves



重量 / Weight 0.22kg

材質 / Material

the airbag is TPU

26kan

海綿雙面復合邊倫布,氣囊為TPU

尺碼 / Size 342*408*194mm

壓力範圍 / Pressure Range

重量 / Weight 0.11kg

left and right : S, M, L,XL,XXL

sponge double-sided composite edge Lunbu

訓

產

品

系

列

20

Pneumatic

Training

Product

Serie

材質 / Material 高分子材料 polymer materials 尺碼 / Size 340*210*70mm **SA-06**

重量 / Weight 1.315kg 分指手腕關節訓練器 Split finger wrist joint trainer



手腕訓練手套

材質 / Material ABS塑料,尼龍軟粘帶 ABS plastic, nylon soft adhesive tape

尺碼 / Size left and right :210*290mm

壓力範圍 / Pressure Range 正壓180kap.負壓80kap P+180kap,P-80kap

Wrist training gloves 重量 / Weight 0.38kg



鏡像手套 Mirror gloves



重量 / Weight 0.06kg

材質 / Material

尺碼 / Size

潛水布料 diving cloth

左手,右手:S,M,L號

left and right : S, M, L

材質 / Material

cotton thread

尺碼 / Size



單管鏡像手套 Single channel mirror gloves

棉線 cotton thread 尺碼 / Size

材質 / Material

左手,右手均碼 left and right, free size



五管鏡像手套 Five way mirror gloves



左手,右手均碼 left and right, free size

棉線



鏡像手柄 Mirror handle

材質 / Material PVC

尺碼 / Size 均碼 free size

重量 / Weight



壓力外控手柄 Pressure external control handle

材質 / Material 塑膠 Plastic

尺碼 / Size

均碼 140*80*38mm

壓力設置範圍 / Pressure setting range Systolic pressure: 0-100KPa (± 10%) Opening pressure: 0-60KPa (± 10%)

波 按

列

Air Wave

Massage

Product

Series

Massage



空氣波按摩系列

Air Wave Massage Series



上肢空氣波壓力按摩儀

Air Wave Pressure Massage Instrument For Upper Limb

Model: JCO

主機+六腔雙手臂按摩套

Host + Six Chamber Massage Coat For Double Upper Limb

▶ 產品特點/Features

- 應用空氣波壓力,反復對肢體進行加壓後再卸壓,促進血液循 環,淋巴循環,緩解手臂肌肉疲勞,减少手臂僵硬等問題。
- 一鍵式操作,簡單、方便,可供個人家庭、醫院、養老院等機構。
- 6個充氣腔室,按摩時間和壓力可調節。
- 配置120W品質電機馬達,力道強勁,使用壽命長。
- 舒適PV材質,透氣柔軟。
- 導氣管三通介面可拆卸不鬆動,天然膠管低溫不變硬。
- 適用人群:年老體弱、運動人士、職業酸痛、久坐群體。
- Apply air wave pressure to repeatedly pressurize the limbs and then relieve the pressure, promote blood circulation, lymphatic circulation, relieve arm muscle fatigue, reduce arm stiffness and
- One-button operation, simple and convenient, can be used for personal homes, hospitals, nursing homes and other
- 6 inflatable chambers, massage time and pressure adjustable.
- Equipped with 120W quality motor, it has strong power and long
- Comfortable PV material, breathable and soft.
- The three-way interface of the air guide tube can be disassembled without loosening, and the natural rubber tube will not harden at low temperature.
- Applicable people: frail elderly, sports people, occupational pains, sedentary groups.



▶ 產品規格/Specifications

主機 / Host

材質 / Material 塑膠 plastic

電源 / Power Supply input AC100-240V, output DC12V, 50Hz

設置 / Setting 時間:0~30分鐘,壓力:30~240mmHg 持續最大壓力:235mmHg

time setting: 0~30 minutes, pressure setting: 30~240mmHg continuous maximum pressure:235mmHg

模式 / Material 單模式循環 diving cloth

運行噪音 / Operating noise

≤65 dB

尺寸 / Size (L*W*H) 重量 / Weight 300*230*120mm 2.2Kg



六腔雙上肢按摩套 **Double Upper Limb Massage Coat-Six Chamber**

防水佈+軟體氣囊 waterproof cloth + soft airbag

產品款式 / Material

雙上肢使用,6充氣腔 for both upper limbs, 6 inflatable cavities

套筒尺:均碼,肩帶長度:130cm sleeve ruler: one size shoulder strap length: 130cm

波

按

系

列

24

下肢空氣波壓力按摩儀

Air Wave Pressure Massage Instrument For Lower Limbs

Model: JCN

主機+六腔雙腿部按摩套

Host + Six Chamber Massage Coat For Double Lower Limb

▶ 產品特點/Features

- 應用空氣波壓力,反復對肢體進行加壓後再卸壓,促進血液循 環,淋巴循環,緩解下肢肌肉疲勞,减少腿部僵硬等問題。
- 一鍵式操作,簡單、方便,可供個人家庭、醫院、養老院等機構。
- 6個充氣腔室,按摩時間和壓力可調節。
- 配置120W品質電機馬達,力道強勁,使用壽命長。
- 舒適PV材質,透氣柔軟。
- 導氣管三通介面可拆卸不鬆動,天然膠管低溫不變硬。
- 適用人群:年老體弱、運動人士、職業酸痛、久坐群體。
- Apply air wave pressure to repeatedly pressurize lower limbs and then relieve the pressure, promote blood circulation, lymphatic circulation, relieve lower limb muscle fatigue and
- One-button operation, simple and convenient, can be used for personal homes, hospitals, nursing homes and other
- 6 inflatable chambers, massage time and pressure adjustable.
- Equipped with 120W quality motor, it has strong power and long
- Comfortable PV material, breathable and soft.
- The three-way interface of the air quide tube can be disassembled without loosening, and the natural rubber tube will not harden at low temperature.
- Applicable people: frail elderly, sports people, occupational pains, sedentary groups.



▶ 產品規格/Specifications

材質 / Material

電源 / Power Supply

time setting: 0~30 minutes, pressure setting: 30~240mmHg continuous maximum pressure:235mmHq

模式 / Material

單模式循環 diving cloth

尺寸 / Size (L*W*H) 300*230*120mm

主機 / Host

塑膠 plastic

input AC100-240V, output DC12V, 50Hz

設置 / Setting

時間:0~30分鐘,壓力:30~240mmHg 持續最大壓力:235mmHg

運行噪音 / Operating noise ≤65 dB

> 重量 / Weight 2.2Kg

腰部空氣波壓力按摩儀

Air Wave Pressure Massage Instrument **For Waist**

Model: JCP

主機+六腔腰部按摩套

Host + Six Chamber Massage Coat For Massage

▶ 產品特點/Features

- 應用空氣波壓力,反復對肢體進行加壓後再卸壓,促進血液循 環,淋巴循環,緩解腰部疲勞酸楚疼痛等問題。
- 一鍵式操作,簡單、方便,可供個人家庭、醫院、養老院等機構。
- 按摩時間和壓力可調節。
- 配置120W品質電機馬達,力道強勁,使用壽命長。
- 舒適PV材質,透氣柔軟。
- 導氣管三通介面可拆卸不鬆動,天然膠管低溫不變硬。
- 適用人群:年老體弱、運動人士、職業酸痛、久坐群體。
- Apply air wave pressure to repeatedly pressurize the limbs and then relieve the pressure, promote blood circulation, lymphatic circulation, relieve waist fatigue, pain and other problems.
- One-button operation, simple and convenient, can be used for personal homes, hospitals, nursing homes and other
- The massage time and pressure can be adjusted.
- Equipped with 120W quality motor, it has strong power and long service life.
- Comfortable PV material, breathable and soft.
- The three-way interface of the air guide tube can be disassembled without loosening, and the natural rubber tube will not harden at low temperature.
- Applicable people: frail elderly, sports people, occupational pains, sedentary groups.



▶ 產品規格/Specifications

主機 / Host

材質 / Material 塑膠 plastic

電源 / Power Supply input AC100-240V, output DC12V, 50Hz

設置 / Setting

時間:0~30分鐘,壓力:30~240mmHg 持續最大壓力:235mmHg

time setting: 0~30 minutes, pressure setting: 30~240mmHg continuous maximum pressure:235mmHa

模式 / Material

單模式循環 diving cloth

運行噪音 / Operating noise ≤65 dB

尺寸 / Size (L*W*H) 300*230*120mm

重量 / Weight 2.2Kg



六腔雙下肢按摩套 Double Lower Limb Massage Coat-Six Chamber

防水佈+軟體氣囊 waterproof cloth + soft airbag

產品款式 / Material 雙上肢使用,6充氣腔 for both lower limbs, 6 inflatable cavities

尺碼 / Size(L*W) 900*320mm



六腔腰部氣動按摩套 Waist Pneumatic Massage Cover - Six Chambers

防水佈+軟體氣囊 waterproof cloth + soft airbag

產品款式 / Material

腰部使用,6充氣腔 for waist use, 6 inflation chambers

尺碼 / Size(L*W) 1200*360mm

波

按

列

Air Wave

上肢空氣波壓力按摩儀

Air Wave Pressure Massage Instrument For Upper Limbs

Model: JCW

主機+四腔單手臂按摩套

Host + Four Chamber Massage Coat For Single Upper Limb

▶ 產品特點/Features

- 應用空氣波壓力,反復對肢體進行加壓後再卸壓,促進血液循 環,淋巴循環,緩解手臂肌肉疲勞,减少手臂僵硬等問題。
- 一鍵式操作,簡單、方便,可供個人家庭、醫院、養老院等機構。
- 4個充氣腔室,按摩時間和壓力可調節。
- 配置120W品質電機馬達,力道強勁,使用壽命長。
- 舒適PV材質,诱氣柔軟。
- 導氣管三通介面可拆卸不鬆動,天然膠管低溫不變硬。
- 適用人群:年老體弱、運動人士、職業酸痛、久坐群體。
- Apply air wave pressure to repeatedly pressurize the limbs and then relieve the pressure, promote blood circulation, lymphatic circulation, relieve arm muscle fatigue, reduce arm stiffness and
- One-button operation, simple and convenient, can be used for personal homes, hospitals, nursing homes and other
- 4 inflatable chambers, massage time and pressure adjustable.
- Equipped with 120W quality motor, it has strong power and long
- Comfortable PV material, breathable and soft.
- The three-way interface of the air guide tube can be disassembled without loosening, and the natural rubber tube will not harden at low temperature.
- Applicable people: frail elderly, sports people, occupational pains, sedentary groups.



▶ 產品規格/Specifications

主機 / Host

材質 / Material 塑膠 plastic

電源 / Power Supply input AC100-240V, output DC12V, 50Hz

設置 / Setting

time setting: 0~30 minutes, pressure setting: 30~240mmHg continuous maximum pressure:235mmHg

模式 / Material

運行噪音 / Operating noise ≤65 dB

尺寸 / Size (L*W*H) 300*230*120mm

重量 / Weight

時間:0~30分鐘,壓力:30~240mmHg 持續最大壓力:235mmHg

單模式循環 diving cloth

2.2Kg

Upper Limb Massage Coat - Four Chamber



防水佈+軟體氣囊 waterproof cloth + soft airbag

產品款式 / Material

單上肢使用, 4充氣腔 Single upper limb use, 4 inflatable cavities

尺碼 / Size(L*W) 900*320mm

下肢空氣波壓力按摩儀

Air Wave Pressure Massage Instrument For Lower Limbs

Model: JCT

主機+四腔單腿部按摩套

Host + Four Chamber Massage Coat For Single Lower Limb

▶ 產品特點/Features

- 應用空氣波壓力,反復對肢體進行加壓後再卸壓,促進促進血液 循環,淋巴循環,緩解下肢肌肉疲勞,减少腿部僵硬等問題。
- 一鍵式操作,簡單、方便,可供個人家庭、醫院、養老院等機構。
- 4個充氣腔室,按摩時間和壓力可調節。
- 配置120W品質電機馬達,力道強勁,使用壽命長。
- 舒適PV材質,透氣柔軟。
- 導氣管三通介面可拆卸不鬆動,天然膠管低溫不變硬。
- 適用人群:年老體弱、運動人士、職業酸痛、久坐群體。
- Apply air wave pressure to repeatedly pressurize lower limbs and then relieve the pressure, promote blood circulation, lymphatic circulation, relieve lower limb muscle fatigue and reduce leg stiffness.
- One-button operation, simple and convenient, can be used for personal homes, hospitals, nursing homes and other
- 4 inflatable chambers, massage time and pressure adjustable.
- Equipped with 120W quality motor, it has strong power and long
- Comfortable PV material, breathable and soft.
- The three-way interface of the air guide tube can be disassembled without loosening, and the natural rubber tube will not harden at low temperature.
- Applicable people: frail elderly, sports people, occupational pains, sedentary groups.



▶ 產品規格/Specifications

主機 / Host

材質 / Material 塑膠 plastic

電源 / Power Supply input AC100-240V, output DC12V, 50Hz

設置 / Setting 時間:0~30分鐘,壓力:30~240mmHg 持續最大壓力:235mmHg


time setting: 0~30 minutes, pressure setting: 30~240mmHg continuous maximum pressure:235mmHa

模式 / Material 單模式循環 diving cloth

運行噪音 / Operating noise ≤65 dB

尺寸 / Size (L*W*H) 300*230*120mm

重量 / Weight 2.2Kg



四腔單下肢按摩套 Upper Limb Massage Coat - Four Chamber

防水佈+軟體氣囊 waterproof cloth + soft airbag

產品款式 / Material

單下肢下使用, 4充氣腔 single lower limb, 4 inflatable cavities

尺碼 / Size(L*W) 900*320mm

波

按

系

列

28

腰部空氣波壓力按摩儀

Air Wave Pressure Massage Instrument For Waist

Model: JCU

主機+四腔腰部按摩套

Host + Four Chamber Waist Massage Coat

▶ 產品特點/Features

- 應用空氣波壓力,反復對肢體進行加壓後再卸壓,促進血液循 環,淋巴循環,緩解腰部疲勞酸楚疼痛等問題。
- 一鍵式操作,簡單、方便,可供個人家庭、醫院、養老院等機構。
- 按摩時間和壓力可調節。
- 配置120W品質電機馬達,力道強勁,使用壽命長。
- 舒適PV材質,透氣柔軟。
- 導氣管三通介面可拆卸不鬆動,天然膠管低溫不變硬。
- 適用人群:年老體弱、運動人士、職業酸痛、久坐群體。
- Apply air wave pressure to repeatedly pressurize the limbs and then relieve the pressure, promote blood circulation, lymphatic circulation, relieve waist fatigue, pain and other problems.
- One-button operation, simple and convenient, can be used for personal homes, hospitals, nursing homes and other
- The massage time and pressure can be adjusted.
- Equipped with 120W quality motor, it has strong power and long
- Comfortable PV material, breathable and soft,
- The three-way interface of the air guide tube can be disassembled without loosening, and the natural rubber tube will not harden at low temperature.
- Applicable people: frail elderly, sports people, occupational pains, sedentary groups.



▶ 產品規格/Specifications

主機 / Host

材質 / Material 塑膠 plastic

電源 / Power Supply

設置 / Setting

time setting: 0~30 minutes, pressure setting: 30~240mmHg continuous maximum pressure:235mmHg

模式 / Material 單模式循環 diving cloth

運行噪音 / Operating noise

尺寸 / Size (L*W*H) 300*230*120mm

重量 / Weight 2.2Kg

input AC100-240V, output DC12V, 50Hz

時間:0~30分鐘,壓力:30~240mmHg 持續最大壓力:235mmHg

≤65 dB

空氣波壓力腿部按摩儀

Air Wave Pressure Leg Massage Instrument

Model: MAR

主機+大腿按摩套+小腿連腳套

Host + Thigh Massage Cover + Calf Foot Cover

▶ 產品特點/Features

- 材質柔軟舒適,氣壓仿真人按摩,釋放腿部壓力。
- 氣動按摩原理,依次充氣放氣迴圈按摩。
- 通過氣囊的收縮擠壓按捏肌肉,放鬆繃緊的肌肉。
- 波動揉捏通過氣壓可按個人需求調整達到按壓舒適力度。
- 9檔力度模式調節,胖瘦皆宜,根據自身的承受能力進行不同檔 位自由選擇。
- 適用人群:需要健康美腿、勞站工人、長期上班缺乏鍛煉、久坐久 站、處於亞健康狀態的人群。
- The material is soft and comfortable, and the air pressure imitates the real person massage to release the pressure on the
- Pneumatic massage principle, successively inflate, deflate and
- Squeeze and pinch the muscles through the contraction of the air bag to relax the tight muscles.
- Fluctuating kneading can be adjusted according to personal needs through air pressure to achieve comfortable pressing
- 9 gear force mode adjustment, suitable for both fat and thin, and different gears can be freely selected according to their
- Applicable crowd: people who need healthy legs, labor station workers, lack of exercise at work for a long time, sedentary and standing for a long time, and are in sub-health status.



▶ 產品規格/Specifications

主機 / Host

材質 / Material 塑膠 plastic

電源 / Power Supply input AC100-240V, output DC12V, 50-60Hz

額定功率 / Rated Power

額定工作時間 / Rated Working Time 15~250 minutes

力度調節 / Strength Adjustment 9檔力度 9 levels of strength

電源線長度 / Power Cord Length 1500mm

重量 / Weight 1.5Kg



四腔腰部氣動按摩套


Waist Pneumatic Massage Cover - Four Chambers

防水佈+軟體氣囊 waterproof cloth + soft airbag

產品款式 / Material

腰部使用,4充氣腔 for waist use, 4 inflation chambers

尺碼 / Size(L*W) 1200*360mm



腿部按摩套 Leg Massage Coat

材質 / Material

防水佈+軟體氣囊 waterproof cloth + soft airbag

產品款式 / Material

小腿和足部,大腿使用 calf, foot, and thigh use

尺寸 / Size 均碼 Average code

空氣波壓力腿部按摩儀

Air Wave Pressure Leg Massage Instrument

Model: MAS

主機+腿部按摩套

Host + Leg Massage Cover

▶ 產品特點/Features

- 空氣波壓縮腿部按摩,可緩解小腿疲勞和疼痛,改善睡眠,增加 血液循環。
- 大腿按摩可以緩解長期坐著和站立引起的血液循環。
- 加熱膝蓋,增加血液循環,驅寒除濕,緩解膝蓋疼痛。
- 智能定時控制,10分鐘自動按摩時間。
- 兩種按摩模式選擇:

順序按摩:大腿和下小腿按摩。

整體按摩:大腿和下小腿同時按摩。

- Air wave compression leg massage can relieve leg fatigue and pain, improve sleep and increase blood circulation.
- Thigh massage can relieve the blood circulation caused by long-term sitting and standing.
- Heat the knee, increase blood circulation, expel cold and dehumidification, and relieve knee pain.
- Intelligent timing control, 10 minute automatic massage time.
- Two massage modes:

Sequential massage: massage thighs and lower legs in turn.

Overall massage: massage the thigh and lower leg at the same time.



▶ 產品規格/Specifications

主機 / Host

材質 / Material 塑膠 Plastic

電源 / Power Supply input AC100-240V, output DC12V, 50-60Hz

額定功率 / Rated Power 熱敷溫度 / Temperature 24W 40~60°C

額定工作時間 / Rated Working Time 10 minutes

力度調節 / Strength Adjustment 3檔力度: 35kPa,45kPa,50kPa 3 levels of strength: 35kPa,45kPa,50kPa

電源線長度 / Power Cord Length 1500mm

尺碼 / Size(L*W*H) 重量 / Weight 215*65*51mm 0.998Kg



腿部按摩套 Leg Massage Coat

> 材質 / Material 蜂窩透氣面料 honeycomb breathable fabric

使用腿圍 / Use Leg Circumference 小腿尺寸:30-48cm,大腿尺寸:40-60cm calf size: 30-48cm, thigh size: 40-60cm

尺碼 / Size(L*W*H) 560*300*400mm